

High-power ultrasound for reprocessing in clinics and practices





Cleaning and chemical disinfection
of medical instruments

Content


SONOREX	
Ultrasonic baths for cleaning and chemical disinfection of medical instruments	1
Recommendations on the application	2
Knowledge of ultrasound.....	3
SONOREX Ultrasonic baths	
Digital or analogue – compact or built-in bath – your choice!	4
SONOREX DIGITEC DT	
Ultrasonic baths with digital operation	5
SONOREX SUPER RK	
Ultrasonic baths with easy-to-operate turning knobs	5
SONOREX ZE	
Ultrasonic bats for built-in.....	6 – 7
SONOREX	
Accessories	8 – 9
SONOMIC	
Ultrasonic baths for MIS instruments and standard instruments	10 – 12
Accessories and Consumables	13
TRISON	
Ultrasonic baths for robotic instruments, MIS instruments and standard instruments	14 – 17
Accessories and Consumables	18
SONOBOARD Sets	
Ultrasonic baths in stainless steel cabinets	19 – 20
Accessories	21
The Foil test.....	21
STAMMOPUR	
Cleaning and disinfection agents	22 – 23
Company portrait.....	24

SONOREX

Ultrasonic baths for cleaning and chemical disinfection of medical instruments



Medical instruments before ...



... and after cleaning with ultrasound

Ultrasound and cavitation – what are they?

Oscillations at frequencies above 18 kHz (18,000 oscillations per second) are termed ultrasound. During the tension phase these oscillations cause the generation of millions of tiny vacuum bubbles in all liquids, which then implode in the compaction phase, causing highly effective pressure surges. This process is called cavitation.

How does ultrasound help in the cleaning of instruments?

Cavitation causes dirt residue and infectious agents in the instruments placed in the liquid to be positively "exploded away". Locations, surfaces, corners and openings that are hard to access can be reached and cleaned by ultrasound ("electronic brushes"). Ultrasound performs the cleaning in a few minutes, and surpasses any manual cleaning.

At the same time it acts gently, as it causes no mechanical damage such as scratching. BANDELIN ultrasonic baths operate usually at the intensive cleaning frequency of 35 kHz. All are equipped with Sweep to provide a homogeneous ultrasound field.

Advantages of instrument treatment in the ultrasonic bath

- High cleaning effect in hard-to-access locations such as drill holes, articulations or joints – with no mechanical damage
- Gentle on instruments
- Rapid instrument circulation
- Reduction of chemical disinfection (duration) by catalytic effect when using suitable preparations (e.g. STAMMOPUR DR 8).
- Economical use of water, chemicals and energy

Recommendations on the application

BANDELIN ultrasonic baths, in combination with the right accessories and preparations made especially for use with ultrasonic baths, permit fast and thorough cleaning and chemical disinfection of medical instruments.

They are used

- as mechanical support for manual cleaning processes
- for removing stubborn soiling before or after mechanical reprocessing
- for cleaning support as an integral part of the mechanical reprocessing procedure.
- for shorter disinfection times while still maintaining intensive cleaning levels

It is important to remember that all cleaning objects must be thoroughly rinsed under running water after use in the ultrasonic bath.



Objects to be cleaned	Cleaning and disinfection agents	Usage notes
Standard instruments (scissors, needle holders, tweezers, forceps, trocars)	STAMMOPUR R Intensive cleaning STAMMOPUR DR 8 VAH certified, for simultaneous cleaning and chemical disinfection	Direct sonication in insert baskets following dry deposits or non-fixing wet deposits Silicone knob mats for placement of sensitive instruments
MIS instruments and accessories, micro-clamps, take-apart tube shaft instruments		
Micro-instruments for neurosurgery/ophthalmology		
Endoscope accessories such as biopsy forceps, snares, valves		
EKG/EEG elektrodes		
Small parts		
Robotic instruments	STAMMOPUR R Intensive cleaning	Direct sonication , connected at the moving device
Stained, encrusted or oxidised instruments	STAMMOPUR GR Basic cleaning	Indirect sonication in an insert beaker or insert tub

Knowledge of ultrasound

Which ultrasonic bath should I select?

The size of the cleaning object will determine the size of the bath and thus the device type. Basket dimensions must be taken into account when selecting a device. To prevent device overload, it is always better to choose a somewhat larger device. This also results in additional space for other uses.

Further important criteria for the decision are the operating controls and the desired design – see following page. For rinsable MIS instruments and complex robotic instruments, ultrasonic baths with additional functions such as rinsing and instrument motion are available, in order to meet the higher cleaning requirements.

Does an ultrasonic bath need a heating?

Devices without heating are preferred for cleaning after dry deposits, as at temperatures above 40 °C there is a risk of protein coagulation. Devices with a heating are used for basic cleaning of instruments, as in such cases, heating of the bath fluid shortens the cleaning time and removes soiling more quickly.

What accessories are necessary?

Cleaning objects must not lie on the bottom of the tank. Baskets and other insert beakers prevent scratching both to the cleaning objects and the bottom of the tank. When cleaning very small or sensitive parts, further accessories may be advisable to facilitate careful placement. For safety reasons, it is recommended to operate ultrasonic baths covered during operation.

What fluids should be used?

STAMMOPUR preparations have been specially developed for use in ultrasonic baths. Water without a detergent will not have a cleaning effect. Do not use household detergents or pure fully desalinated (DI) water. For work with acids, a plastic insert tub must be used. Never use inflammable or explosive fluids directly in the oscillating tank!

How can ultrasonic baths be tested?

The effectiveness of ultrasonic baths depends on the intensity and distribution of the process-typical



cavitation in the oscillating tank. The foil test (in accordance with DIN 58341:2020-07; IEC TR 60886) is a simple procedure for demonstrating the intensity and distribution of cavitation in an ultrasonic bath. In this test, an aluminium foil placed in the tank is perforated / destroyed to a certain degree by cavitation, depending on the duration of sonication. To achieve reproducible foil test results, it is important to provide similar testing conditions in each case. Suitable setups for performing foil tests are available as accessories for the ultrasonic baths.

More information about foil test see of page 23.

If you want to know more ...

... visit our website or our YouTube channel with a lot of helpful tutorials!

Or contact us directly... we are always pleased to provide advice, so call us at +49 (0)30 76880-212.



SONOREX Ultrasonic baths

Digital or analogue – compact or built-in bath – **your choice!**



	DIGITEC DT ... /M	SUPER RK ...	ZE ... DT	ZE ...
Capacity [l]	0.9 – 46.0	0.9 – 28.0	29.0 – 46.0	29.0 – 46.0
Ultrasonic frequency [kHz]	35	35	35	35
Pulse function	✓	✓	✓	✓
Sweep	✓	✓	✓	✓
Rapid degassing DEGAS	✓	–	✓	–
Additional ultrasound from the side	–	–	ZE 1032 DT / 1059 DT	ZE 1032 / 1059
Time setting [min]	1, 2, 3, 4, 5, 10, 15, 30, ∞	1 – 15, ∞	1, 2, 3, 4, 5, 10, 15, 30, ∞	1 – 15, ∞
Program storage	1 program	–	1 program	–
Safety shut-down	after 12 hours	–	after 12 hours	–
Heating thermostatically adjustable [°C]	version "H" : 20 – 80	version "H" 30 – 80 RK 31 H: 65 fixed	–	–
Setting accuracy of bath temperature [K]	± 2,5	± 5	–	–
Excess temperature signal	✓	–	✓	–
Inclined tank bottom for complete emptying	DT 1058 M	–	✓	✓
Filling level mark for safe dosage	✓	✓	✓	✓
Outlet	one-piece welded from DT 106	one-piece welded from RK 106	bead 1½" (drain set G 1½ optional)	
Thickness tank material [mm]	0.8	0.8	2.0	2.0
Degree of protection	IP 33	IP 32	–	–
Mains supply: 230 V~ (± 10 %), 50/60 Hz 115 V~ (± 10 %), 50/60 Hz	✓ ✓	✓ ✓	✓ ✓	✓ ✓
Medical device	class I	class I	class I	class I

SONOREX DIGITEC DT

Ultrasonic baths with digital operation

SONOREX SUPER RK

Ultrasonic baths with
easy-to-operate
turning knobs



Type	Internal tank dimensions l x w x d [mm]	Capacity [l]	Code No.	External dimensions l x w x d [mm]	Ultrasonic peak power [W]	Ultrasonic nominal power [W]	Heating power [W]	Outlet ball valve
DT 31 DT 31 H RK 31 RK 31 H	190 x 85 x 60	0.9	3200 3220 329 044	205 x 100 x 180	160	40	- 70 - 70	-
DT 100 DT 100 H RK 100 RK 100 H	240 x 140 x 100	3.0	3210 3230 301 312	260 x 160 x 250	320	80	- 140 - 140	-
DT 106 RK 106	Dia. 240 x 130	5.6	3270 306	Dia. 265 x 270	480	120	- -	G ¼
DT 156 RK 156	500 x 140 x 100	6.0	3275 305	530 x 165 x 245	640	160	- -	G ¼
DT 255 DT 255 H RK 255 RK 255 H	300 x 150 x 150	5.5	3215 3240 3066 316	325 x 175 x 295	640	160	- 280 - 280	G ¼
DT 514 DT 514 H RK 514 RK 514 H	325 x 300 x 150	13.5	3250 3211 277 207	355 x 325 x 305	860	215	- 600 - 600	G ½
DT 1028 DT 1028 H RK 1028 RK 1028 H	500 x 300 x 200	28.0	3255 3231 322 324	535 x 325 x 400	1200	300	- 1300 - 1300	G ½
DT 1058 M	600 x 400 x 200/220+	46.0	304120	670 x 470 x 400	2400	600	-	G ¾

*corresponds to 4 times output + inclined tank bottom

SONOREX ZE



Ultrasonic baths for built-in

Advantages of built-in baths

- Hygienic, unobstructed work surfaces thanks to under-table mounting
- Inclined tank bottom for easier emptying
- Hygienic maintenance due to rounded tank corners
- Operating control on the front side
- Ultrasonic generators may be installed optional right or left in the base cupboard
- Suitable for 1/1 DIN baskets as of model ZE 1031 and ISO baskets as of ZE 1058
- Built-in bath with ultrasound and rinsing tank without ultrasound – an option to expand your worktop

Built-in baths with bottom and side sonication

The foil test figures below show that ultrasonic baths with bottom and side sonication generate a more homogeneous sound field than baths with bottom sonication alone. This means a more gentle and uniform cleaning, an important consideration for highly sensitive instruments.



Mounting examples

- Optimum sound distribution and reduction of acoustic shadows as a result of additional side sonication
- Electronically induced movements of the sound field by means of TwinSonic technology reduce local peaks of impact
- No additional oscillation necessary for the instrument basket, and no further space is needed in the working area
- The latest generator technology with Sweep
- Existing built-in baths with bottom sonication are easy to replace, thanks to an identical tank edge design

Type	Internal tank dimensions l × w × d [mm]]	Capacity [l]	Code No.	External dimensions l × w × d [mm]	Ultrasonic peak power [W]	Ultrasonic nominal power [W]	Outlet
------	--	-----------------	----------	--	------------------------------	---------------------------------	--------

... with sonication from bottom

ZE 1031 ZE 1031 DT	510 × 300 × 200/220 ⁺	29.0	3060 3217	570 × 360 × 270/290 ⁺	1200	300	bead 1½"
ZE 1058 ZE 1058 DT	600 × 400 × 200/220 ⁺	46.0	3050 3234	660 × 460 × 270/290 ⁺	2400	600	bead 1½"

... with sonication from bottom and side

ZE 1032 ZE 1032 DT	510 × 300 × 200/220 ⁺	29.0	3075 3223	570 × 404 × 270/290 ⁺	1760	440	bead 1½"
ZE 1059 ZE 1059 DT	600 × 400 × 200/220 ⁺	46.0	3085 3248	660 × 504 × 270/290 ⁺	2400	600	bead 1½"

Rinsing tanks without ultrasound, for mounting into worktops

suitable for bath	Type	Code No.	Internal tank dimensions l × w × d [mm]	External dimensions l × w × d [mm]	Description
ZE 1031/1032 / ... DT	SW 31 Z	3048	510 × 300 × 200/220 ⁺	570 × 360 × 205/225 ⁺	with bead 1½", without drain set
ZE 1058/1059 DT	SW 58 Z	3049	600 × 400 × 200/220*	660 × 460 × 205/225 ⁺	with bead 1½", without drain set
<i>optional: drain set G 1½</i>		3166			<i>drain set with turning knob and stainless steel stopper</i>

Digital control unit with temperature display









suitable for bath	Type	Code No.	Description
ZE 1031 bis ZE 1059	ST 30 DT	309803	The ST 30 DT digital control unit has an integrated temperature display and offers the user added safety to prevent protein coagulation. If the bath fluid heats up to > 40 °C, a red warning LED will also flash.

*corresponds to 4 times output + inclined tank bottom, mounting from below

Dimension without ultrasound generator, external dimensions of ultrasound generators 360 × 310 × 142 mm (l × w × h)

SONOREX Accessories

Appropriate accessories facilitate ultrasonic application and simultaneously protect the oscillating tank and instruments.
Objects to be cleaned or beakers must not be placed onto the tank bottom!

Accessories	Material	Function	Images (selection)	RK 31 / H DT 31 / H	RK 100 / H DT 100 / H
Lid	stainless steel	protects the bath fluid from contaminants from the outside condensation water runs back into the tank lid D...T especially for inset baskets with drop handles	 D 514 Type Code No.	D 08 218	D 100 / D 3 T 3003 / 114
Hinged lid	stainless steel	protects the bath fluid from contaminants from the outside condensation water runs back into the tank hinged lid D...G for built-in units	 D 1031 G Type Code No.	–	–
Insert basket	stainless steel	to use for the instruments to be cleaned	 K 14 Type I x w x d [mm] Code No.	K 08 170 x 65 x 50 209	K 3 C 200 x 110 x 40 3025
Inset basket	stainless steel	to use for the instruments to be cleaned. baskets with hinged handles in standard formats for instrument cleaning without basket change basket holder KT is necessary	 K 29 EM Type I x w x d [mm] Code No.	–	–
Basket holder	stainless steel	support for inset baskets or DIN 1/1 and 1/2 sieve trays KT...Z is equipped with handles	 KT 57 Type Code No.	–	–
Insert tub	plastic, with lid	especially for basic instrument cleaning with STAMMOPUR GR don't use at temperatures higher 60 °C	 KW 3 Type I x w x d [mm] Code No.	–	KW 3 195 x 115 x 88 715
Knob mat	silicone	for contact-free placement of highly-sensitive instruments, especially micro-instruments, during cleaning prevents damage to instruments; permeable for ultrasound	 SM 14 Type I x w [mm] Code No.	–	SM 3 170 x 97 093
Frame for foil test	stainless steel	The frame is used for foil test, which is as simple method for displaying the intensity and distribution of the cavitation in an ultrasonic bath.	 FT 1 Type Code No.	FT 1 3190	FT 4 3074

RK 106 DT 106	RK 156 DT 156	RK 255 / H DT 255 / H	RK 514 / H DT 514 / H	RK 1028 / H DT 1028 / H	DT 1058 M	ZE 1031 / DT ZE 1032 / DT	ZE 1058 / DT ZE 1059 / DT
D 6 346	D 156 3004	D 255 / D 5 T 3007 / 3054	D 514 / D 14 T 3010 / 3062	D 1028 / D 28 T 3011 / 3063	D 1058 M 7526	D 30 7522	D 57 7520
–	–	–	–	–	–	D 1031 G 3229	D 1058 G 3232
K 6 dia. 215 x 50 356	K 6 L 460 x 100 x 50 202	K 5 C 260 x 110 x 40 3027	K 14 275 x 245 x 50 354	K 28 455 x 245 x 50 358	–	–	–
–	–	–	–	K 29 EM 470 x 240 x 45 688	K 29 EM 470 x 240 x 45 688	K 29 EM 470 x 240 x 45 688	K 29 EM 470 x 240 x 45 688
–	–	–	–	KT 30 7517	KT 57 7504	KT 30 / KT 30 Z 7517 / 7507	KT 57 / KT 57 Z 7504 / 3078
–	–	KW 5 254 x 96 x 130 240	KW 14 280 x 215 x 145 613	KW 28-0 437 x 230 x 155 717	–	–	–
–	SM 6 426 x 97 110	SM 5 213 x 97 101	SM 14 235 x 245 118	SM 29 470 x 245 178	SM 29 470 x 245 178	SM 29 470 x 245 178	SM 29 470 x 245 178
FT 4 3074	FT 6 3222	FT 4 3074	FT 14 3084	FT 40 3094	FT 37 3674	FT 36 3673	FT 37 3674

SONOMIC

Ultrasonic baths for MIS instruments and standard instruments

Three patents in one device!



Mounting example SONOMIC MC 1001 E

The reliable internal cleaning of MIS instruments and rinsable parts of other instruments ensures their continued use.

SONOMIC has been specially developed for these instruments and combines the effects of

damage-free ultrasonic cavitation,
the effective suction rinsing and
individual testing of instruments

in one device.

The integrated flow-control monitoring for each connected instrument guarantees reliable cleaning results and prevents instrument malfunction.

Advantages at a glance:

- Safety as a result of patented individual instrument rinsing and testing
- Patented suction rinsing principle
- Patented universal adapter for instrument connection without change of seal
- Temperature and filling level monitoring
- Reproducible program sequences
- Versatility:
 - Can be used for standard instruments too
- Documentation by means of protocol print-outs
- Available as compact bath or built-in bath


Individual instrument examination rather than general testing

If different MIS instruments are rinsed at the same time, the rinsing result for the individual instruments cannot be checked.

With SONOMIC this problem is solved by means of the innovative channel selector (patent EP 19 20 797). Only one instrument at a time is released for rinsing, thus permitting individual flow-through monitoring. The minimum flow-through rate is 2 ml/s at 0.8 bar vacuum, otherwise the instrument will be identified as "non rinsable" and so indicated on the touch screen. The determination, classification and clear indication of successful rinsability for each instrument ensures a higher safety level during reprocessing.


SONOMIC suction rinsing compared to pressure rinsing from other suppliers

Generally, the majority of the soiling is concentrated at the distal end of MIS instruments. In comparable devices from other manufacturers, all MIS instruments are rinsed under pressure from the distal end. As a result, contamination is forced through the whole lumen of the instruments, thus presenting an increased risk of undesirable deposits, especially in constricted areas near the handles and in other difficult-to-reach areas of the instruments. The direct entry of dissolved contamination into the bath fluid is a further negative effect. The suction rinsing function (patent DE 20 2006 020 523) exerted by the SONOMIC at the distal end of the instruments removes soiling against the direction of penetration, while fresh cleaning solution takes its place. This avoids unnecessary contamination of the rear lumen parts of the instruments. The removed contamination moves through the adapter into the exchangeable filter, rather than back into the bath fluid.



Connection of instruments to the universal adapter without change of seal

In the SONOMIC, twelve rinsable MIS instruments with diameters from 3 to 10 mm can each be connected to one of the identical adapters without having to change the adapter seal. The seal has an innovative torsion principle (patent EP 19 20 727) that guarantees complete fluid-tightness against the outer shaft of the instruments. The highly elastic sealing material has been ultrasonically tested and is resistant to the preparations STAMMOPUR R and STAMMOPUR DR 8. A maintenance-related exchange of the seals is necessary after more than 500 load cycles. Because of the user-friendly adapter design the exchange of the seals can be applied without any tools.



Torsion principle of the adapter seal

Filling level and temperature monitoring

The SONOMIC monitors the correct bath fluid level with an integrated filling level sensor. In case of non-compliance, it will not be possible to start the sonication, and the user will receive an error message on the touch screen.

Before each cycle, a temperature sensor tests whether the bath temperature is within the permitted range. If the bath fluid temperature is too low (< 18 °C), the heating automatically switches on. To prevent protein coagulation, a warning message appears at temperatures of about 40 °C.

Safety as a result of reproducible program sequences

The SONOMIC operating program contains a self-test and provides the user with clear instructions for all necessary work stages. For instance, an adapter test is performed prior to each load, an indispensable measure for reliable identification of non-penetrable instruments.

Documentation by means of protocol print-outs

For quality verification, SONOMIC provides several interfaces. When required, protocols with the following data can be printed out: cleaning mode, bath temperature, result of rinsing examination, etc.

SONOMIC MC 1001

Ready-to-use set:

- Ultrasonic bath MC 1001
- Basket K 1001 MC
- 12 adapters with seal and hose ADS 1000
- 12 adapter seals AD 1000
- Adapter testing strips APB 1000
- 30 filter cartridges EF 1001
- Frame for foil test FT 38



SONOMIC MC 1001 E for built-in

Consisting of:

- Oscillating tank TE 1001 E
- Ultrasound generator GT 1001 E
- Control unit ST 1001 E
- Basket K 1001 MC
- 12 adapters with seal and hose ADS 1000
- 12 adapter seals AD 1000
- Adapter testing strips APB 1000
- 30 filter cartridges EF 1001
- Frame for foil test FT 38



Internal tank dimensions l x w x d [mm]	650 x 400 x 160/170 ⁺	650 x 410 x 160/170 ⁺
Capacity [l]	42.5	43.5
Operating volume [l]	27.0	27.5
Code No.	3315	3345
External dimensions l x w x d [mm]	860 x 490 x 325	oscillating tank: 855 x 475 x 250 ultrasound generator (with rinsing module): 485 x 155 x 360
Ultrasonic peak power* [W]	2400	2400
Ultrasonic nominal power [W]	600	600
Ultrasonic frequency [kHz]	40	40
Control: touchscreen 96 x 61 mm	✓	✓
Heating, program-controlled [W]	400	400
Pulse function	✓	✓
Sweep	✓	✓
Temperature monitoring	✓	✓
Thickness tank material [mm]	2.0	2.0
Filling level mark	✓	✓
Level sensor	✓	✓
Outlet	ball valve ¾", thread feed pipe G ¾", at the rear side	G 1½ drain set with turning knob and stainless steel stopper
Mounting into the worktop	–	from below
Interfaces	USB-B, RS-232, LPT	USB-B, RS-232, LPT
Medical device	class I	class I





*corresponds to 4 times output + inclined tank bottom

SONOMIC

Accessories and Consumables

Accessories		MC 1001	MC 1001 E
Lid Code No.		D 1000 MC 3312	
Hinged lid Code No.		D 1001 G 3310	D 1001 GE 3326
Inset basket l x w x d [mm] Code No.		K 1001 MC 520 x 340 x 50 3324	
Knob mat l x w [mm] Code No.		SM 1000 MC 245 x 340 3313	
Frame for foil test Code No.		FT 38 550 x 470 3672	

Consumables

Filter cartridges Code No.		EF 1001 à 30 pcs. 3365	EF 1001 à 100 pcs. 3366	
Adapter seals Code No.		AD 1000 à 12 pcs. 3353	AD 1000 à 24 pcs. 3354	AD 1000 à 36 pcs. 3355
Adapters with seals and hose Code No.		ADS 100 à 1 pc. 3350	ADS 100 à 12 pcs. 3351	
Adapter testing strips Code No.		APB 1000 à 2 pcs. 3358		

TRISON

Ultrasonic baths for robotic instruments, MIS instruments and standard instruments

The new standard for pre-cleaning of robotic instruments



Mounting example TRISON Xi

An effective cleaning process for robotic instruments verifiably takes place when constant movement allows optimum access by the ultrasound to the working tools and cables of the instrument.

TRISON is a new modular ultrasonic bath for intensive pre-cleaning of high-grade medical instruments, particularly robotic surgery instruments.

For the first time, and unique in the world, it combines

Ultrasound
Individual rinsing and
Moving of instruments

for optimum cleaning results in complex robotic instruments.

Advantages at a glance:

- Improved cleaning through a combination of ultrasound, rinsing and moving for robotic instruments
- Designed for robotic instruments
- Reliability as a result of individual instrument rinsing and examination
- Simple instrument connection
- For robotic instruments, but also for rinsable MIS and standard instruments
- Temperature monitoring with warning function
- Protocol function
- Ergonomic, hygienic controls
- Flexible, space-saving system for the workplace thanks to different mounting options

Moving improves the cleaning effect!

The heart of the TRISON is the innovative moving device TRISON Twist for up to four robotic instruments. Special sprung actuators engage with the bottom of the instruments and move them during sonication. The surgical working tools are rotated and manipulated at their tips for all-over ultrasonic effect. Integrated friction clutches prevent damage to filigree cables and drive rollers.



Instrument tip in movement

The cleaning process is also supported by the permanent internal rinsing of the instruments. To allow this, the instruments are connected with hoses to the rinsing cycle of the TRISON Base control unit.

The removed soiling is guided directly into the exchangeable filter, rather than back into the bath fluid. Various series of experiments with actually contaminated instruments and with standardised test contamination confirm the effectiveness of the new TRISON cleaning concept.

Ultrasonic bath specifically for robotic instruments

The TRISON ultrasonic oscillating bath has been specially dimensioned for extremely long robotic instruments. Thanks to the high-performance ultrasonic oscillating systems on the bottom and sides, drive adapters, instrument shafts and working tools at their tips are optimally reached by ultrasound and possible acoustic shadows are avoided. The TRISON ultrasound generator is equipped with Sweep in order to minimise standing waves and to guarantee a homogeneous ultrasonic intensity distribution. In a cleaning program designed in consultation with the manufacturer, robotic instruments are first soaked for approx. 30 min. in order to partially dissolve or break down organic residues. In the subsequent alternating suction and pressure rinsing process with ultrasound support, the control housing and shaft of each instrument is rinsed and tested for flow-through. In this way, even

stubborn contamination is reliably stripped off and removed.

Individual instrument rinsing and examination for even greater reliability


During sonication, each instrument is individually released for rinsing and checked for flow-through (patent DE 20 2006 020 523). Non-rinsable instruments are reliably identified and displayed on the touch screen at the end of the process. The determination, classification and clear indication of successful rinsability for each instrument ensure a higher level of safety for reprocessing.

Simple instrument connection

The TRISON Twist allows the fixation of up to four robotic instruments by means of a simple push-on mechanism – no additional basket is required. Connection to the rinsing cycle is made with one double Luer connector per instrument.

The special TRISON Rack basket allows up to eight MIS instruments to be connected for rinsing, using the reliable SONOMIC adapter.

Connection to the TRISON Base control unit is performed quickly and without risk of mix-up by using two multi-hose connectors.



Robotic instruments before ...



... and after cleaning



TRISON Twist for robotic instruments



TRISON Rack for MIS instruments



Inset basket for standard instruments

Versatility: Three cleaning options in one device

TRISON has been specially developed for robotic instruments. With the use of suitable accessories, however, MIS instruments and standard instruments can also be effectively cleaned. In order to best meet the various cleaning requirements, a cleaning program is available for each instrument type.

Temperature monitoring with warning function

In the modes for cleaning of MIC instruments and robotic instruments, an included sensor tests whether the temperature of the bath fluid is within the permitted range. For the cleaning of standard instruments in the inset basket, a separate temperature sensor is required.


Depending on the used cleaning agent, the operator can define temperature limits to ensure an effective cleaning process.

Logging of the reprocessing cycles through an ethernet or USB interface

For quality verification, TRISON provides two interfaces. A USB interface allows for simple and flexible data transfer using a USB stick. The ethernet interface allows link-up to central sterile services management software or similar. Logging includes inter alia: cleaning mode, bath temperature, result of rinsing examination.

Ergonomically designed user interface, hygienic touch screen operation

When designing TRISON, special value was assigned to ergonomics and fitness for use in the intended work environment within a hospital's central sterile services department. The result was device control via an easy-to-clean touch screen. The user-friendly operator software includes many self-explanatory images without lengthy operating texts.



Flexible installation options

The swivelling control unit TRISON Base can be flexibly mounted/installed on the right or left side of the ultrasonic oscillating tank. This allows TRISON to be adapted to the different installation conditions according to the space available.

TRISON 4000

Consisting of:
















- Oscillating tank TRISON TE 3000 with drain set ①
- Ultrasound generator TRISON GT 3000 M-C ②
- Control unit TRISON Base TB 4000 ③
- Mains supply switch NW 3000 ④
- Moving device TRISON Twist TT 4000 Xi ⑤ or TT 4000 Si ⑥
- Pivot mounted arm TRISON Lift TL 4000 ⑦
- 30 filter cartridges EF 1001
- Frame for foil test FT 42



	TRISON 4000 Xi		TRISON 4000 Si	
Internal tank dimensions, l x w x d [mm]	770 x 420 x 165/190 ⁺			
Capacity [l]	60.0			
Operating volume [l]	35.0			
Code No. (version)	7885 (right)	7985 (left)	7884 (right)	7984 (left)
External dimensions, l x w x d [mm]	oscillating tank: 900 x 480 x 245/275 ⁺ / ultrasound generator: 360 x 310 x 142 control unit: 370 x 190 x 380 / mains supply switch: 220 x 60 x 145 pivot mounted arm: 240 x 95 x 350			
	moving device Xi: 345 x 160 x 175		moving device Si: 405 x 205 x 190	
Ultrasonic peak power* [W]	3040			
Ultrasonic nominal power [W]	760			
Ultrasonic frequency [kHz]	38			
Sonication from side	✓			
Pulse function	✓			
Sweep	✓			
Temperature monitoring	✓			
Thickness tank material [mm]	2.0			
Filling level mark	✓			
Outlet	G 1½ drain set with turning knob and stainless steel stopper			
Mounting into the worktop	from below			
Mains supply: 230 V~ (± 10%), 50/60 Hz	✓			
Current consumption** [A]	3.5			
Interfaces	USB, Ethernet			
Inlet pressure [bar]	5 ... 9, ISO 8573-1 (7:4:4)			
Medical device	class I			

*corresponds to 4 times ultrasonic nominal output ** in case of 230 V~[± 10%] 50/60 Hz ⁺inclined tank bottom

TRISON Accessories and Consumables

Accessories		Type	Code No.		External dimensions l × w × d [mm]	Function and Compatibility
 Moving device TRISON Twist Xi / Si	TT 4000 Xi	right 7821	left 7921	345 × 160 × 175 405 × 205 × 190	for Xi- or Si-robotic instruments available either as right or left version; to use with TRISON Lift	
	TT 4000 Si	7820	7920			
 Pivot mounted arm TRISON Lift	TL 4000	7930		240 × 95 × 350	for TRISON Twist 4000	
 TRISON Rack	TR 3001	right 7631	left 7731	640 × 405 × 150	basket with connections for up to eight MIS instruments; available either as right or left version	
 Inset basket	K 29 EM	688				
 Basket holder	KT 3000 Z	7761		–	support of the inset basket	
 Lid	D 4000 A	right 7955	left 7956	–	universal for all TRISON applications; for TRISON Twist only in lowered position; available either as right or left version	
 Spacer	TX 4000 Xi	7763		–	for TRISON Twist Xi, PU à 2 pcs.	
 Frame for foil test	FT 42	3224		700 × 440	The frame is used for foil test, which is as simple method for displaying the intensity and distribution of the cavitation in an ultrasonic bath.	
 Temperature sensor	TM 4000	7741		–	for measuring the temperature of the bath during the cleaning of standard instruments	
Consumables						
Filter cartridges Code No.		EF 1001 à 30 pcs. 3365			EF 1001 à 100 pcs. 3366	
Adapter seals Code No.		AD 1000 à 8 pcs. 3361			AD 1000 à 24 pcs. 3354	
Adapter with seal Code No.		ADT 1000 à 1 pc. 7770			ADT 1000 à 8 pcs. 3359	
Adapter testing strips Code No.		APB 3000 à 2 pcs. 7771				
Hose set with coupling for TRISON Twist Code No.		for Xi: SLS 4000 TT à 1 pc. 3362			for Si: SLS 3000 TT à 1 pc. 3363	
Hose set with coupling for TRISON Rack, without adapters Code No.		SLS 3000 TR 3364				

SONOBOARD Sets

Ultrasonic baths in stainless steel cabinets

The practical supplement to your sink unit facility, or for individual use!



SONOBOARD TRISON

For selected ultrasonic baths, BANDELIN offers ready-to-use sets consisting of an ultrasonic bath and a practical stainless steel cabinet. The double-walled stainless steel cabinets are equipped with overlapping fronts and all-round rubber seals on the doors and panels. Their flexible positioning (thanks to lockable casters), ergonomic working height and additional storage space make them a high quality and irreplaceable item of clinic equipment.

SONOBOARD has a high resilience to scratches and impacts, and is extremely resistant against chemicals. The smooth stainless steel surfaces reduce the accumulation of germs and bacteria to a minimum, and meet the most rigorous hygiene requirements.

■ Start-up and operation is fast and simple!

SONOBOARD is delivered as a ready-to-use set; only the utility services need to be connected. Three variations are available, designed for different instrument types.

Features SONOBOARD STANDARD

- Operation on the front side
- Digital control unit with temperature monitoring
- Suitable for 1/1 DIN and ISO baskets

Features SONOBOARD MIC

- Safety as a result of patented individual instrument rinsing and testing
- Patented suction rinsing principle
- Patented universal adapter for instrument connection without change of seal
- Temperature and filling level monitoring
- Reproducible program sequences

Features SONOBOARD TRISON

- Improved cleaning through a combination of ultrasound, rinsing and movement for robotic instruments
- Designed for robotic instruments
- Reliability as a result of individual instrument rinsing and testing
- Simple instrument connection
- Temperature monitoring (Robotic and MIS mode)
- Protocol function
- Available as left and right version

SONOBOARD STANDARD for standard instruments

Ready-to-use set:

- Ultrasonic bath ZE 1058 with accessories (see page 8 – 9)
- Functional cabinet FS 900 S



SONOBOARD MIC for MIS instruments and standard instruments

Ready-to-use set:

- Ultrasonic bath SONOMIC MC 1001 E with accessories (see page 13)
- Functional cabinet FS 1200 ML



SONOBOARD TRISON for robotic instruments, MIS and standard instruments

Ready-to-use set:

- Ultrasonic bath TRISON 4000 Xi with accessories (see page 18)
- Functional cabinet FS 1200 T



	SONOBOARD STANDARD	SONOBOARD MIC	SONOBOARD TRISON	
Internal tank dimensions l x w x d [mm]	600 x 400 x 200 / 220 ⁺	650 x 410 x 160 / 170 ⁺	770 x 420 x 165/190 ⁺	
Capacity [l]	46.0	43.5	60.0	
Operating volume [l]	32.0	27.5	35.0	
Code No. (version)	3452	3454	3457 (right)	3456 (left)
External dimensions incl. rolls, l x w x h [mm]	900 x 700 x 930	1200 x 700 x 930	1200 x 700 x 930	
Ultrasonic peak power* [W]	2400	2400	3040	
Ultrasonic nominal power [W]	600	600	760	
Ultrasonic frequency [kHz]	35	40	38	
Pulse function	✓	✓	✓	
Sweep	✓	✓	✓	
Time setting [min]	1, 2, 3, 4, 5, 10, 15, 30, ∞	menu controlled	menu controlled	
Temperature monitoring	✓	✓	✓	
Thickness tank material [mm]	2.0	2.0	2.0	
Filling level mark	✓	✓	✓	
Level sensor	–	✓	–	
Outlet	G 1½ drain set with turning knob, stainless steel stopper and siphon G 1½ with hose			
Mains supply: 230 V~ (± 10 %), 50/60 Hz	✓	✓	✓	
115 V~ (± 10 %), 50/60 Hz	✓	✓	–	
Current consumption** [A]	2.7	3.0	3.5	
Interfaces	–	USB-B, RS-232, LPT	USB-A, Ethernet	

*corresponds to 4 times output **in case of 230 V ~ (± 10 %) 50/60 Hz ⁺ inclined tank bottom

SONOBOARD

Accessories


	SONOBOARD STANDARD		SONOBOARD MIC		SONOBOARD TRISON	
for standard instruments Code No.	Basket holder KT 57 Z 3078	Insert basket K 29 EM 688	–	–	Basket holder KT 3000 Z 7761	Insert basket K 29 EM 688
for MIS instruments Code No. (version)	–	–	Insert basket K 1001 MC 3324	–	TRISON Rack TR 3001 7631 (right) 7731 (left)	–
for robotic instruments Code No. (version)	–	–	–	–	TRISON Twist TT 4000 Xi 7821 (right) 7921 (left)	TRISON Twist TT 4000 Si 7820 (right) 7920 (left)
Lid Code No. (version)	D 57 7520	–	D 1000 MC 3312	–	D 4000 A 7955 (right) 7956 (left)	–
Hinged lid Code No.	D 1058 G 3232	–	D 1001 GE 3326	–	–	–
Frame for foil test Code No.	FT 37 3674	–	FT 38 3672	–	FT 42 3224	–

The Foil test – Function testing of an ultrasonic bath

A foil test¹ is recommended for testing ultrasonic baths. This should be conducted upon initial startup and at regular intervals thereafter (e.g. every 3 months). The frequency of testing is the responsibility of the user. The foil test is a simple procedure for demonstrating the intensity and distribution of cavitation in an ultrasonic bath. It involves stretching aluminium foil over a foil testing frame, which will be perforated or destroyed to a certain degree by cavitation, depending on sonication time.

For purposes of reproducibility, it is important that the test conditions remain constant:

- Filling the oscillation tank to the filling level mark
- Temperature of the sonication fluid
- Degassing time
- Positioning of frame
- Foil type (brand, thickness)
- Sonication time
- Type and concentration of ultrasonic agent



https://bandelin.com/foil_test/

Foils can be archived in a suitable way (scanning, photos, etc.) This allows the foils to be compared at any time. The perforated areas of all foils should have approx. the same dimensions and distribution – the results are never identical.

A process validation, e.g. for the treatment of medical products, can only be achieved by conducting regular foil tests.

To execute the foil test, different foil test frames can be ordered from the manufacturer BANDELIN. The foil test frames are suitable for a wide range of tank dimensions.

Aluminium household foil is also required to conduct the test and is not included in the delivery.

STAMMOPUR

Cleaning and disinfection agents

For optimum cleaning results in the ultrasonic bath, specially formulated detergents and disinfection agents are required alongside ultrasound performance, temperature and time. BANDELIN offers a balanced range of cleaning agents and disinfectants from DR. H. STAMM GmbH.

With their cavitation-conductive properties, these preparations support the cleaning process while protecting the materials.

The preparations are biologically degradable in accordance with the regulations of the Detergents Directive. In treating the instruments, it is important to rinse them thoroughly after using the ultrasonic bath.



Important: Some common cleaning and disinfection agents from other manufacturers may contain components that attack the ultrasonic oscillating tank and could lead to breakdown due to pitting.

Product information and EC safety data sheets are available as pdf downloads from safetydatasheets.bandelin.com

For the best cleaning results, the indicated concentration of the preparation must be observed. To facilitate dosing, we provide a dosing table available that is only suitable for Bandelin equipment. The dosage table is available online: dosingtable.bandelin.com

Unit	Operating volume (l/min)	1%	2%	3%	5%	10%
DR 3014, DR 3014	0,8	500 ml	1000 ml	1500 ml	2500 ml	5000 ml
DR 3014, DR 3014	1,2	750 ml	1500 ml	2250 ml	3750 ml	7500 ml
DR 3014, DR 3014	2,0	1250 ml	2500 ml	3750 ml	6250 ml	12500 ml
DR 3014, DR 3014	3,0	1875 ml	3750 ml	5625 ml	9375 ml	18750 ml
DR 3014, DR 3014	4,0	2500 ml	5000 ml	7500 ml	12500 ml	25000 ml
DR 3014, DR 3014	6,0	3750 ml	7500 ml	11250 ml	18750 ml	37500 ml
DR 3014, DR 3014	8,0	5000 ml	10000 ml	15000 ml	25000 ml	50000 ml
DR 3014, DR 3014	10,0	6250 ml	12500 ml	18750 ml	31250 ml	62500 ml
DR 3014, DR 3014	12,0	7500 ml	15000 ml	22500 ml	37500 ml	75000 ml
DR 3014, DR 3014	15,0	9375 ml	18750 ml	28125 ml	46875 ml	93750 ml
DR 3014, DR 3014	20,0	12500 ml	25000 ml	37500 ml	62500 ml	125000 ml
DR 3014, DR 3014	30,0	18750 ml	37500 ml	56250 ml	93750 ml	187500 ml
DR 3014, DR 3014	40,0	25000 ml	50000 ml	75000 ml	125000 ml	250000 ml
DR 3014, DR 3014	60,0	37500 ml	75000 ml	112500 ml	187500 ml	375000 ml
DR 3014, DR 3014	80,0	50000 ml	100000 ml	150000 ml	250000 ml	500000 ml
DR 3014, DR 3014	100,0	62500 ml	125000 ml	187500 ml	312500 ml	625000 ml
DR 3014, DR 3014	150,0	93750 ml	187500 ml	281250 ml	468750 ml	937500 ml
DR 3014, DR 3014	200,0	125000 ml	250000 ml	375000 ml	625000 ml	1250000 ml
DR 3014, DR 3014	300,0	187500 ml	375000 ml	562500 ml	937500 ml	1875000 ml
DR 3014, DR 3014	400,0	250000 ml	500000 ml	750000 ml	1250000 ml	2500000 ml
DR 3014, DR 3014	600,0	375000 ml	750000 ml	1125000 ml	1875000 ml	3750000 ml
DR 3014, DR 3014	800,0	500000 ml	1000000 ml	1500000 ml	2500000 ml	5000000 ml
DR 3014, DR 3014	1000,0	625000 ml	1250000 ml	1875000 ml	3125000 ml	6250000 ml
DR 3014, DR 3014	1500,0	937500 ml	1875000 ml	2812500 ml	4687500 ml	9375000 ml
DR 3014, DR 3014	2000,0	1250000 ml	2500000 ml	3750000 ml	6250000 ml	12500000 ml
DR 3014, DR 3014	3000,0	1875000 ml	3750000 ml	5625000 ml	9375000 ml	18750000 ml
DR 3014, DR 3014	4000,0	2500000 ml	5000000 ml	7500000 ml	12500000 ml	25000000 ml
DR 3014, DR 3014	6000,0	3750000 ml	7500000 ml	11250000 ml	18750000 ml	37500000 ml
DR 3014, DR 3014	8000,0	5000000 ml	10000000 ml	15000000 ml	25000000 ml	50000000 ml
DR 3014, DR 3014	10000,0	6250000 ml	12500000 ml	18750000 ml	31250000 ml	62500000 ml

Dosing table

Preparation	Description	Application with ultrasound Concentration, Duration	Litres	Code No.
STAMMOPUR R CE - Concentrate -	Intensive cleaner of medical instruments in the ultrasonic bath. High cleaning efficiency, even for instruments heavily contaminated with incrustations of blood and secretions. Anticorrosive, very high material compatibility, applicable for all materials. Also applicable as contact liquid in the ultrasonic bath – e.g. for recommended basic cleaning of spotted and ugly looking instruments with STAMMOPUR GR. Without phosphates, aldehydes and chlorine. Main active agents: tensides, mildly alkaline, pH 9.6 at 1 %.	2 %, 2 – 10 min	2 5 10	934 989 1029
STAMMOPUR DR 8* – VAH-certified – CE 0483 - Concentrate -	Manual, chemical disinfection and disinfecting, non-fixing intensive cleaning of general, surgical, invasive and non-invasive medical instruments and accessories in an ultrasonic bath as well as in an immersion bath process. High blood dissolution, for instruments heavily contaminated with incrustations of blood and secretions. Short irradiation time. Solution applicable under strain for 3 sequent days. Very high material compatibility. Concentrate. Non-odiferous. Without aldehydes, chlorine, phenols. Bactericidal, yeasticidal, limited virucidal, additionally active against H5N1, SV40, Adeno. Mildly alkaline pH 9.4 at 1 %. Labelling in accordance with CLP. Signal word: Danger, GHS05-GHS07-GHS08-GHS09	2 %, 5 min SV40 with high protein burden: 2 %, 10 min Adeno with high protein burden: 3 %, 15 min application without ultrasound: 1 %, 60 min: bactericidal, levurocidal; 2% – 30 min or 3% – 15 min: bactericidal, levurocidal, limited virucidal incl. H5N1 and additionally against SV40	2 5 10	972 974 1028
STAMMOPUR GR* CE - Concentrate -	Basic cleaning of spotted, encrusted and ugly looking instruments in the ultrasonic bath. Removes tarnish, metal oxides, rust, spotting, burned-in residues after sterilisation and mineral residues e.g. lime. Caution with damaged chroming and nickel-plated parts. Not for light metals, tin and zinc. Not to be used for routine cleaning. Application at 50 – 60 °C only in a insert tub. Main active agents: phosphoric acid, tensides, pH 1.9 at 1 %. Only to be used for basic cleaning.	5 %, 2 – 10 min	2 5 10	938 969 1031

*Transport regulations for 'Dangerous goods' have to be observed. No transport by airfreight.

Dosing aids	usable for	Code No.
Pump ①	5-l-jerrycan	268
Pump ①	10-l-jerrycan	2660
Measuring beaker ②	100 ml	294



BANDELIN Ultrasound since 1955

Company portrait

We are a family-owned company located in Berlin and meanwhile run in the third generation, specialised in development, manufacturing and sales of ultrasonic devices, the corresponding accessories and application-specific cleaning agents and disinfectants.

A wide vertical range of manufacture, modern production lines and a motivated staff guarantee a high quality of the products. Our devices contribute to the success of our customers in the laboratory, medical, dental, pharmaceutical, industrial, craft as well as service.

As early as 1955, our company began developing and manufacturing high-performance ultrasonic devices. The constant expansion of the product range and a sharp rise in sales led to an expansion of the production area in 1985. In 1992, ultrasonic homogenisers and controllable, power-constant ultrasonic generators were introduced to the market.

The period from 1996 to 2004 was characterised by the development and production of innovative ultrasonic baths and immersible transducers as well as tube reactors for industrial applications.

In the following years, BANDELIN's product range was expanded by new laboratory ultrasonic devices. After the introduction of the ultrasonic bath for simultaneous cleaning and rinsing of MIS instruments, a further development was launched in 2016 for robotic instruments.

Today, the reputation of our brands SONOREX, SONOPULS, SONOMIC and TRISON stand for the high quality awareness of our employees and is equated in expert circles with ultrasound.

The most important product groups include:

- SONOREX – ultrasonic baths and reactors
- SONOPULS – ultrasonic homogenisers
- SONOMIC – ultrasonic baths for rinsable MIS and standard instruments
- TRISON – ultrasonic baths for robotic-, rinsable MIS and standard instruments
- TICKOPUR – cleaning agents
- STAMMOPUR – cleaning agents and disinfectants

We are innovation leaders in the development of ultrasonic devices and new areas of application. In the past we have registered 79 patents / utility models as well as 68 trade brands. Our participation in various committees in the development of new standards and guidelines serve to ensure the highest standards for ultrasonic applications.

As the only complete supplier of ultrasonic devices, accessories, cleaning agents and disinfectants with approvals and certifications according to ISO 9001 and ISO 13485, BANDELIN is the market leader.

Over one million units have already been delivered to our customers.



More information about our company you will find here: bandelin.com/prospekte/Company_history_GB.pdf

Made in Germany

BANDELIN electronic
GmbH & Co. KG
Heinrichstraße 3 – 4
12207 Berlin
DEUTSCHLAND
☎ +49 30 76880-0
📠 +49 30 7734699
info@bandelin.com

Certified in accordance with
ISO 9001 and ISO 13485



We advice you personally!
Feel free to consult our experts:

+49 30 76880-212

www.bandelin.com



6912-002 EN/2021-09

Subject to technical alterations without notice.

Trademarks: da Vinci, da Vinci Si, da Vinci Xi and
Intuitive Surgical are registered trademarks
of Intuitive Surgical Inc.

All images are provided as examples and are not true to size.
Decorative elements are not included in the scope of delivery.

The general terms and conditions apply.